

**KCMC Biotechnology  
Laboratory, Microbiology**

**STANDARD  
OPERATING  
PROCEDURE**

**Effective Date**  
**09 January 07**

**SOP-Number**  
**MIC.041**

**Page 1 of 7**

**Date**  
**09 January 2007**

**Title: EXAMINATION OF BLOOD SMEARS FOR MALARIA PARASITES**

**Author/Date: Anne Morrissey, December 2006**

**Approvals/Date:**

\_\_\_\_\_  
**John A. Crump, MB, ChB, Microbiology Laboratory Director**

\_\_\_\_\_  
**Anne Morrissey, MS, MT(ASCP)SM, Microbiology Lab Supervisor**

**This SOP has been read and understood by:**

| <b>Name</b>                         | <b>Date</b> |
|-------------------------------------|-------------|
| <b>Signatures on original copy.</b> |             |
| <b>1.</b>                           |             |
| <b>2.</b>                           |             |
| <b>3.</b>                           |             |
| <b>4.</b>                           |             |
| <b>5.</b>                           |             |
| <b>6.</b>                           |             |
| <b>7.</b>                           |             |
| <b>8.</b>                           |             |
| <b>9.</b>                           |             |
| <b>10.</b>                          |             |

| <b>Annual Review</b> |             |
|----------------------|-------------|
| <b>By</b>            | <b>Date</b> |
|                      |             |
|                      |             |
|                      |             |
|                      |             |
|                      |             |

**Title: EXAMINATION OF BLOOD SMEARS FOR MALARIA PARASITES**

**Document History:**

| <b>Version Number</b> | <b>Reason for Changes</b> | <b>Date</b> |
|-----------------------|---------------------------|-------------|
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |
|                       |                           |             |

**Copies distributed to:**

| <b>Name</b> | <b>Date</b> |
|-------------|-------------|
|             |             |
|             |             |
|             |             |

**Title: EXAMINATION OF BLOOD SMEARS FOR MALARIA PARASITES**

**PURPOSE**

Procedures for the examination of thick and thin Giemsa stained blood smears for the presence and identification of *Plasmodium* spp. (*P. falciparum*, *P. vivax*, *P. ovale*, *P. malariae*).

**PRINCIPLE**

Species specific identification of the *Plasmodium* spp. is clinically important for several reasons: *P. falciparum* has a high mortality rate and more likely to be drug resistant, eradication of *P. vivax* and *P. ovale* infections requires supplemental therapy. The thick blood film with lysed red blood cells (RBCs) is approximately 30 times more sensitive than the thin film for detecting the presence of malaria parasites (detecting about 20 parasites/ $\mu$ l). The thin film with intact RBCs is examined for speciation of malarial parasites allowing evaluation of the position of organisms within the cells, the size and shape of infected erythrocytes and presence of stippling within the cells.

**SPECIMEN COLLECTION AND HANDLING**

Refer to SOP MIC.040 PREPARATION AND STAINING OF SMEARS FOR BLOOD PARASITE EXAMINATION.

**BLOOD PARASITE SMEAR RESULT LOG**

Enter the following information on the Result Log:

Study ID#

Patient Name

KCMC Medical Record Number

Specimen Collection Date

**QUALITY CONTROL**

**Positive Malaria Slide:**

**Frequency:** Each new batch of stain prior to use and each day of patient testing.

**Control:** Thin smears prepared from blood from patient with high parasitemia (one parasite per 2-3 hpf), fixed in methanol and stored at -80°C.

**Acceptable Results:**

1. Background should be clear and free from debris.
2. Nuclei of WBCs stain deep purple, cytoplasm pale purplish-blue.
3. RBCs stain pale pink.
4. Malaria parasite – blue rings, red chromatin dot.

**Corrective Action for Unacceptable Results:**

1. Restain new slide under supervision and reexamine.

**Title: EXAMINATION OF BLOOD SMEARS FOR MALARIA PARASITES**

2. If still unsatisfactory, prepare fresh stain and repeat QC.

**Documentation:**

1. Record QC results on GIEMSA STAIN QC SHEET.
2. Document all corrective actions for unacceptable results on QC Deviation form and submit to supervisor for review.
3. Review QC results monthly.

**Slide Review:**

1. As part of ongoing training and competency assessment of staff, all slides will be reviewed within 24 hours of initial reading with emphasis on suspected positive smears.
2. Discordant results between initial examiner and reviewer will be reviewed by the with technician and documented on Blood Parasite Smear Result Log.

**EXAMINATION OF THICK SMEARS**

1. Examine entire thick smear using low or high power objectives (10X, 25X) to detect large blood parasites such as microfilariae. Refer to SOP MIC.042 IDENTIFICATION OF NON-MALARIA BLOOD PARASITES if any are suspected. If other blood parasite is suspected have the slide reviewed immediately by Lab Supervisor.
2. Examine 100 fields using oil immersion objective (100X) for malaria parasite morphology. Select areas that are well stained and free of precipitate.
3. If no parasites are observed in 100 fields, continue examination until 300 fields have been observed.

**EXAMINATION OF THIN SMEARS**

1. Thin smears may be available earlier than thick smears which require longer drying time so can be examined first. If this is the case proceed as follows:
  - a. Screen at low magnification (10X, 20X) for the presence of large blood parasites and proceed as instructed for thick smears.
  - b. Examine 100 fields using oil immersion objective (100X). If no parasites are observed in 100 fields, continue examination until 300 fields have been observed
2. If thick smear has been examined and is negative, examination of the thin smear is not necessary.
3. If parasites are observed, quantify the number of asexual (trophozoites, schizonts) and sexual (gametocyte) forms (See Section – **QUANTIFYING PARASITES**).

**Title: EXAMINATION OF BLOOD SMEARS FOR MALARIA PARASITES**

**IDENTIFICATION OF MALARIA PARASITES**

The identification of the various *Plasmodium* species should be approached in an orderly fashion based on the observation of the following key morphologies:

- Ring forms and other asexual stages found in circulating in blood
- Appearance of the RBCs – enlarged, irregular or stippled
- Size and shape of developing and mature trophozoites
- Number of merozoites per schizont
- Morphology of gametocytes and internal pigment.

Refer to Table 1 and reference textbooks for examples of differentiating morphologies of *Plasmodium* spp. and infected erythrocytes.

**Table 1. Differentiating features of *Plasmodium* spp.**

| Species                      | Stages found in peripheral blood  | Appearance of RBC          |                 | Appearance of Parasite  |  |  |
|------------------------------|---|----------------------------|-----------------|---|--|--|
|                              |   | Size                       | Stippling       | Cytoplasm   | Pigment in mature schizonts, gametocytes | No. of merozoites/schizont                       |
| <i>Plasmodium falciparum</i> | Rings (mature)<br>Gametocytes   | Normal                     | None            | Double dot in rings common<br>Rings often small and delicate<br>>1 ring /RBC                    | Black                                    | 6-32 (schizonts rarely seen in peripheral blood) |
| <i>Plasmodium vivax</i>      | Trophozoites (all stages)<br>Schizonts<br>Gametocytes   | Enlarged                   | Schuffners dots | Ring forms in all stages<br>Ameboid trophs<br>Light blue<br>Large, circular gametocytes         | Golden brown                             | 12-24  |
| <i>Plasmodium ovale</i>      | Trophozoites (all stages)<br>Schizonts<br>Gametocytes   | Enlarged, irregular shapes | Schuffners dots | Rounded compact trophs<br>Dark to medium blue, usually dense                                    | Dark brown                               | 6-14   |
| <i>Plasmodium malariae</i>   | Trophozoites (all stages)<br>Schizont<br>Gametocytes<br>Few rings usually seen –ring stage short. | Normal                     | None            | Rounded compact trophs<br>Dark blue dense cytoplasm<br>Trophozoite band forms seen occasionally | Dark brown coarse                        | 6-12   |

**Title: EXAMINATION OF BLOOD SMEARS FOR MALARIA PARASITES**

**QUANTIFYING MALARIA PARASITES**

The number and type of morphological forms observed (asexual –ring forms, trophozoites, schizonts vs. sexual – gametocytes) is an important diagnostic tool for appropriate treatment. They may indicate possible drug resistance and/or determine prognosis if serial specimens are examined after treatment is initiated.

Several different schemes for quantifying malaria parasites exist. In order to maintain consistency the procedure for quantifying parasites from Tanzanian patients will be the same as is conventionally used.

1. On the thick smear, using the 100x objective count the number of asexual forms (trophozoites, schizonts) seen per 200 WBCs.
2. If any sexual forms (gametocytes) are observed during the examination of the fields above, continue to examine fields containing 300 more WBCs and count the number per 500 WBCs.

**PROCEDURAL NOTES**

1. Care must be taken to avoid identifying superimposed platelets as malarial parasites. Superimposed platelets routinely exhibit a "HALO" whereas malarial parasites do not. It is also helpful to compare the suspected parasite to platelets that normally appear between the erythrocytes.

**REFERENCE RANGE**

Normal - No parasites found

**RESULT REPORTING-MALARIA PARASITES**

If a malarial parasite is seen but identification to species level is difficult, results should be reported: "*Plasmodium* sp. seen". If it is possible to determine that it is not *P. falciparum*, report "*Plasmodium* sp., not *P. falciparum*".

1. Enter test results BLOOD PARASITE SMEAR REPORT LOG.
2. Fill out a MICROBIOLOGY/IMMUNOLOGY RESULT FORM and ISAAC MISC TEST RESULTS DATA FORM and enter test results.
3. A positive slide for a blood parasite is considered a **Critical Result** and must be reported as soon as possible to physician taking care of the patient.
4. Document receipt of test information in Critical Results Log book and have physician sign the log book to acknowledge receipt of information.
5. Place negative patient reports in ward Lab Results box.

**Title: EXAMINATION OF BLOOD SMEARS FOR MALARIA PARASITES**

**REFERENCES**

MalariaGEN: Guidelines for Diagnostic Procedures for Blood Samples, Centers for Disease Control, Atlanta, GA.

Diagnostic Procedures for Blood Specimens for, Centers for Disease Control, Atlanta, GA, July 2003.

Rogers, W.O. 2003. *Plasmodium* and *Babesia*. In: Manual of Clinical Microbiology. 8<sup>th</sup> Ed. ASM Press, Washington, DC.

Centers for Disease Control, Division of Parasitic Diseases website: [www.dpd.cdc.gov/DPDx](http://www.dpd.cdc.gov/DPDx)

APPENDIX A - Blood Parasite Smear Results Log

APPENDIX B – Giemsa Stain QC Record

APPENDIX C – ISAAC Microbiology Test Report Form