

**KCMC Biotechnology  
Laboratory, Microbiology**

**STANDARD  
OPERATING  
PROCEDURE**

**Effective Date**

**SOP-Number  
MIC.039**

**Page 1 of 8**

**Date**

**Title: GIEMSA STAIN FOR BLOOD PARASITE SMEARS**

**Author/Date: Anne Morrissey, December 2007**

**Approvals/Date:**

**John A. Crump, MB, ChB, Microbiology Laboratory Director**

**Anne Morrissey, MS, MT(ASCP)SM, Microbiology Lab Supervisor**

**This SOP has been read and understood by:**

<b>Name</b>	<b>Date</b>
<b>Signatures on original copy.</b>	
<b>1.</b>	
<b>2.</b>	
<b>3.</b>	
<b>4.</b>	
<b>5.</b>	
<b>6.</b>	
<b>7.</b>	
<b>8.</b>	
<b>9.</b>	
<b>10.</b>	

<b>Annual Review</b>	
<b>By</b>	<b>Date</b>

**Title: GIEMSA STAIN FOR BLOOD PARASITE SMEARS**

**Document History:**

<b>Version Number</b>	<b>Reason for Changes</b>	<b>Date</b>

**Copies distributed to:**

<b>Name</b>	<b>Date</b>

**Title: GIEMSA STAIN FOR BLOOD PARASITE SMEARS**

**PURPOSE**

Procedures for preparation and staining of thick and thin blood smears for examination for blood parasites.

**PRINCIPLE**

Whole blood is collected in a vacutainer tube containing EDTA anticoagulant. A large drop of blood is spread in a small circle on a clean, grease-free microscope slide and allowed to air dry (thick smear). A second slide is prepared with a small drop of blood and spread thinly over the slide surface (thin smear). The thin smear will be fixed with methanol prior to staining to preserve the red cell morphology. The thick smear is not fixed and the RBCs lyse during the staining process. Both slides are stained with Giemsa and examined for malaria and other blood parasites. The thick blood film is about 30 times more sensitive than the thin film (detecting about 20 parasites/ $\mu$ l) and is most suitable for the rapid detection of blood parasites, particularly when they are few. The thin film with intact RBCs is examined for speciation of malarial parasites if this is not apparent on the thick smear.

**SCOPE**

This Standard Operating procedure applies to all personnel who prepare and stain slides for blood parasites and are trained and competent in performing these procedures.

**STANDARD PRECAUTIONS**

Standard precautions should be observed when obtaining and handling patient samples, preparing and handling slides. Wear gloves to prevent exposure to bloodborne pathogens.

**SPECIMEN COLLECTION AND HANDLING**

**Specimen Type:** Whole blood collected in EDTA vacutainer tube. Refer top MIC.012 SPECIMEN COLLECTION ISAAC 002 and 003 STUDIES for blood collection procedure.

**Specimen Handling**

Refer to SOP MIC.011 SPECIMEN RECEIPT AND HANDLING for minimal volumes, rejection criteria, pre and post analysis handling and storage.

**Title: GIEMSA STAIN FOR BLOOD PARASITE SMEARS**

### **SMEAR PREPARATION**

#### **Materials:**


2 glass slides cleaned with 70% alcohol to remove grease and film.  
Third slide for spreading thin smear ("spreader slide")  
Pasteur pipette  
Pencil  
Slide box

#### **PROCEDURE**

1. Label 2 slides with the following:  
    Study ID#  
    Date of collection
2. Clean the surface of each slide with 70% alcohol to remove dust and grease.
3. Mix blood sample by gently inverting 5 times.
4. Thick smear
  - a. Using a pasteur pipette dispense one large drop (6 ul) of blood in the center of one slide and spread to a circle with approximately a 1.5 cm diameter.
  - b. Allow to air dry in a horizontal position.
  - c. Place slide in covered slide box.
5. Thin smear
  - a. On the second slide place a small drop (2 ul) of blood in the area adjacent to the frosted end.
  - b. Hold the "spreader slide" at a 35-40° angle up to the drop, and allow the drop to spread along the contact line of the 2 slides.
  - c. Quickly push the spreader slide toward the other (unfrosted) end of the slide.
  - d. Make sure the spread film has good, feathered edge. Repeat if necessary.
  - e. Allow the smear to air-dry with the slide in a horizontal position.
  - f. Place slides in covered slide box.

**Title: GIEMSA STAIN FOR BLOOD PARASITE SMEARS**

**Diagram of smear appearance:**



**Thin Smear**

**Thick smear**

**PREPARATION OF GIEMSA STAIN STOCK SOLUTION**  
**(Alternatively Giemsa Stock stain may be purchased commercially)**

**Materials:**

Giemsa stain powder  
Glycerol  
Methanol (absolute)  
Measuring cylinder  
Funnel & filter paper  
Balance for weighing ingredients  
Mixer  
Brown bottle

**Procedure (For 1 liter of stock stain solution):**

1. Measure 500 ml of glycerol in a graduated cylinder and pour into dark bottle.
2. Add 500 ml absolute methyl alcohol.
3. Add 7.6 grams Giemsa powder.
4. Allow the contents to mix overnight using a mixing machine.
5. Once thoroughly mixed, filter mixture using the filter paper and funnel.
6. Store the filtered solution in brown bottle in a cupboard to avoid exposure to direct sunlight.
7. Clearly label the bottle with the date of preparation, date of expiration 6 months from date of preparation and the initials of the preparer.

**Title: GIEMSA STAIN FOR BLOOD PARASITE SMEARS**

### **PREPARATION OF BUFFER SOLUTION pH 7.2**

**Materials:**

Sodium phosphate, monobasic dehydrate ( $\text{NaH}_2\text{PO}_4$ )  
Potassium phosphate monobasic ( $\text{KH}_2\text{PO}_4$ )  
Distilled water  
Measuring cylinder  
pH meter  
Glass bottle  
Weighing balance

**Preparation:**

1. Weigh 1.09 g  $\text{NaH}_2\text{PO}_4$
2. Weigh 0.7 g  $\text{KH}_2\text{PO}_4$ .
3. Dissolve both in 1000ml distilled water
4. Measure pH. using pH meter.
5. If necessary, add small quantities of 0.1N HCl or 0.1 N NaOH to bring the PH of the buffer to pH 7.2
6. Put solution in glass bottle.
7. Autoclave at  $121^\circ\text{C}$  for 15 min.
8. Label with date of preparation, expiration date 1 year from date of preparation and initials of preparer.
9. Store at room temperature away from sunlight.

### **PREPARATION OF WORKING STAIN SOLUTION**

*Prepare fresh each day of use.*

**Materials:**

Giemsa stock solution  
Buffer solution. pH. 7.2  
Measuring cylinder

**Preparation:**

1. Prepare a 3% solution of the Stock Giemsa stain in phosphate buffer (e.g. For 100 ml of working stain solution add 3 ml of stock Giemsa stain to 97 ml of buffer.)
2. Discard unused solution at the end of each day.

**Title: GIEMSA STAIN FOR BLOOD PARASITE SMEARS**

### **GIEMSA STAIN PROCEDURE**

1. Place THIN smear on staining rack and fix with methanol for 1 minute or immerse in Copelin jar for 1 minute. Allow to air dry. ***DO NOT FIX THICK SMEAR WITH METHANOL.***
2. Place both thick and thin smears on a staining rack and flood slide with working stain solution or fill Copelin jar with the working stain solution and add slides.
3. Stain for 30-45 min.
4. Rinse both smears gently with buffer to remove stain or dip in Copelin jar containing buffer 3-4 times.
5. Flood THICK smear with buffer or immerse in Copelin jar with buffer for 5 minutes.
6. Allow to air dry.

### **QUALITY CONTROL**

**Frequency:** Each new batch of stain.

**Control:** Thin smears prepared from blood from patient with high parasitemia (one parasite per 2-3 hpf), fixed in methanol and stored at -80°C.

#### **Acceptable Results:**

1. Background should be clear and free from debris.
2. Nuclei of WBCs stain deep purple, cytoplasm pale purplish-blue.
3. RBCs stain pale pink.
4. Malaria parasite – blue rings, red chromatin dot.

#### **Corrective Action for Unacceptable Results:**

1. Restain new slide under supervision and reexamine.
2. If still unsatisfactory, prepare fresh stain and repeat QC.

#### **Documentation:**

1. Record QC results on GIEMSA STAIN QC SHEET.
2. Document all corrective actions for unacceptable results on QC Deviation form and submit to supervisor for review.
3. Review QC results monthly.

### **REFERENCES**

MalariaGEN: Guidelines for Diagnostic Procedures for Blood Samples, Centers for Disease Control, Atlanta, GA.

**Title: GIEMSA STAIN FOR BLOOD PARASITE SMEARS**

Diagnostic Procedures for Blood Specimens for, Centers for Disease Control, Atlanta, GA, July 2003.

APPENDIX A - Giemsa Stain QC sheet