

KCMC Biotechnology Laboratory	STANDARD OPERATING PROCEDURE	Effective Date 02-10-2006	SOP-Number FLOW007.01
		Page 1 of 9	02-10-2006
Title: Specimen Accessioning			
SOP References:		Supersedes:	

Author/Date: Moses Sichangi, M.S
Approvals/Date: <div style="margin-left: 100px;"> <hr/> Chris Drakeley, Ph.D. KCMC Biotechnology Laboratory Director </div> <div style="margin-left: 100px;"> <hr/> John A. Bartlett, M.D. ISAAC Program Director </div> <div style="margin-left: 100px;"> <hr/> Guido Ferrari, M.D. ISAAC Pathogenesis Program Co-Director </div>

This SOP has been read and understood by:

Name	Date
1.	
2.	
3.	
4.	
5.	
6.	
7.	
8.	
9.	
10.	
11.	
12.	

KCMC Biotechnology Laboratory	STANDARD OPERATING PROCEDURE	Effective Date 02-10-2006	SOP-Number FLOW007.01
		Page 2 of 9	02-10-2006
Title: Specimen Accessioning			
SOP References:		Supersedes:	

Annual Review		
Reviewed by:	Review Date	Signature

Document History:

Version Number	Reason for Changes	Date
FLOW007.01	Initial	16-08-2006

KCMC Biotechnology Laboratory	STANDARD OPERATING PROCEDURE	Effective Date 02-10-2006	SOP-Number FLOW007.01
		Page 4 of 9	02-10-2006
Title: Specimen Accessioning			
SOP References:		Supersedes:	

1. Method Summary

Before the laboratory can accept a sample, certain criteria must be followed to assure the validity of the specimen.

All specimens accepted by the laboratory must be examined to ascertain that they meet the proper criteria for processing and data entry. In addition, correct processing and data entry is extremely important to assure that the results given out are correct. And finally, there must be verification that the results are for the intended patient.

2. Scope

This SOP applies to all staff at the KCMC Biotechnology laboratory who use the FACSCalibur. Only authorized staff may use the FACSCalibur.

3. Safety Precautions

- 3.1. Follow Standard biosafety Precautions for handling of blood specimens (See Biohazard SOP [SFT.01]).
- 3.2. Any forms that are sent in contact with the specimens (i.e. wrapped around specimen container with rubber band) should be discarded after the proper information is obtained. They are not to be filed.

4. Sample Verification

- 4.1. All specimens are to be accompanied by a completely filled sample collection form and tracking form (see appendix 1 and 2 respectively). The specimen should contain the following information:
 - 4.1.1. Patient ID, Date and time of sample collection, study ID and Visit number (time-point in study).
- 4.2. Forms should be checked to make sure all appropriate information is provided and that it matches the information on the specimen container. Forms should contain the following information:
 - 4.2.1. Patient ID number
 - 4.2.2. Study number (if applicable)

KCMC Biotechnology Laboratory	STANDARD OPERATING PROCEDURE	Effective Date 02-10-2006	SOP-Number FLOW007.01
		Page 5 of 9	02-10-2006
Title: Specimen Accessioning			
SOP References:		Supersedes:	
<div>4.2.3. Time-point in study</div> <div>4.2.4. Clinician requesting for the test</div> <div>4.2.5. Date of sample collection</div> <div>4.2.6. Time sample was drawn</div> <div>4.2.7. The name of the phlebotomist</div> <div>4.2.8. Name and signature of procuring individual (research nurse or research assistants).</div> <div>4.2.9. Sample type and assay type requested for.</div> <div>5. Criteria and procedure for Rejection of Unacceptable Specimens</div> <div>5.1 Specimens may be rejected for the following reasons:<div>5.1.1. Quantity not sufficient (QNS).</div><div>5.1.2. Clotted specimen</div><div>5.1.3 Lack of procurement information as listed on sample collection form (appendix 1).</div><div>5.1.4 Hemolysed blood specimen.</div><div>5.1.5 Sample more than 24 hours old.</div><div>5.1.6 Specimens transported on ice/ice packs.</div><div>5.1.7 Incorrectly labeled specimen (not matching the information on sample forms).</div><div>5.1.8 Unlabeled specimen</div><div>5.1.9 Specimen drawn into tube containing anticoagulant different from protocol specification (i.e EDTA).</div></div> <div>5.2. Specimen Rejection Procedure</div>			
Page 5 of 9			

KCMC Biotechnology Laboratory	STANDARD OPERATING PROCEDURE	Effective Date 02-10-2006	SOP-Number FLOW007.01
		Page 6 of 9	02-10-2006
Title: Specimen Accessioning			
SOP References:		Supersedes:	

5.2.1. Unlabeled or mislabeled specimen should be handled as follows:

- Do not accession.
- Contact procuring individual. If in doubt about labeling information, the sample should be discarded, documented, and patient rescheduled for another sample to be drawn.

5.2.2 Fill out specimen rejection form (see appendix 3).

5.2.3 Notify appropriate clinician and procuring individual.

5.2.4 File specimen rejection form in specimen rejection file.

6. Computer Data Accession Procedure

6.1. The accession book is the data entry point for all samples. No sample should be processed without being assigned an accession number and a daily number. All information received with a specimen should be logged into the database. The forms should be taken to the data management unit of KCMC. The information should be handled according to the data management unit policy # 1.04 for data entry.

Reference:
N/A

Page 6 of 9

KCMC Biotechnology Laboratory	STANDARD OPERATING PROCEDURE	Effective Date 02-10-2006	SOP-Number FLOW007.01
		Page 7 of 9	02-10-2006
Title: Specimen Accessioning			
SOP References:		Supersedes:	

Appendix 1: Sample Collection Form

FORM 1: SAMPLE COLLECTION

Study: _____ **Visit:** _____

1. Date: _____/_____/_____ 2. Sample ID: _____
dd/ mm /yyyy

3. Age: _____ 4. Date of Birth: _____/_____/_____
dd /mm /yyyy

5. Gender: _____ (m/f)

6. Individual consent obtained: _____

7. If yes, HIV antibody test results: () Negative () Positive

8. Volume of blood sample drawn: _____mL (In EDTA anticoagulant)

9. Date and Time of sample collection: _____:_____
hh : mm

10. Name of personnel who obtained sample: _____

Immunology Laboratory

11. Sample received by: _____

12. Date and Time received: _____/_____/_____. _____:_____
dd/ mm/ yyyy hh : mm

Comments:

**KCMC Biotechnology
Laboratory**

**STANDARD
OPERATING
PROCEDURE**

Effective Date
02-10-2006

SOP-Number
FLOW007.01

Page 8 of 9

02-10-2006

Title: Specimen Accessioning

SOP References:

Supersedes:

Appendix 2: Sample Tracking Form

FORM 2: SAMPLE TRACKING

DATE: ____/____/____ STUDY: ____
dd / mm / yyyy

SAMPLE ID	NO. OF VIALS	CONFIRM RECEIPT (√)

Form completed by: Name _____ Sign _____

Date and Time samples released: ____/____/____. ____:____
dd / mm /yyyy hh : mm

Name of Research nurse/assistant: _____ Sign: _____

Immunology Laboratory

Sample received by: _____

Date and time received: ____/____/____. ____:____
dd / mm /yyyy hh :mm

Comments: _____

