

KCMC Biotechnology Laboratory, Microbiology	STANDARD OPERATING PROCEDURE	Effective Date 18 April 07	SOP-Number MIC.043
		Page 1 of 5	Date 18 April 2007
Title: OXIDASE TEST			
Author/Date: Anne Morrissey, 2006			
Approvals/Date: <div style="text-align: center;"> John A. Crump, MB, ChB, Microbiology Laboratory Director </div> <div style="text-align: center;"> Anne Morrissey, MS, MT(ASCP)SM, Microbiology Lab Supervisor </div>			

This SOP has been read and understood by:

Name	Date
Signatures on original copy.	
1.	
2.	
3.	
4.	
5.	
6.	
7.	
8.	
9.	
10.	

Annual Review	
By	Date

KCMC Biotechnology Laboratory, Microbiology	STANDARD OPERATING PROCEDURE	Effective Date 18 April 07	SOP-Number MIC.043
		Page 2 of 5	Date 18 April 2007
Title: OXIDASE TEST			

Document History:

Version Number	Reason for Changes	Date

Copies distributed to:

Name	Date

Title: OXIDASE TEST

PURPOSE

For the detection of cytochrome c oxidase activity in a variety of Gram negative organisms.

PRINCIPLE

A drop of oxidase reagent (N,N,N,N-tetramethyl-p-phenylenediamine dihydrochloride) is placed on a small piece of filter paper. A few colonies of test organism from a non-selective medium (e.g. Blood or Chocolate Agar) are rubbed into the reagent on the filter paper. If cytochrome oxidase is produced by the organism it oxidizes cytochrome c which in turn oxidizes the reagent to form a purple compound.

SCOPE

This Standard Operating Procedure applies to the testing for cytochrome oxidase by technical staff in the microbiology laboratory who have been trained and are competent in the use and performance of this test.

SPECIMEN COLLECTION – NA

EQUIPMENT/MATERIALS

Oxidase reagent (BD) Store at 15-30 C
Filter paper (small piece)
Glass slide
Inoculating loop or wooden applicator stick.
Fresh culture (24-48 hrs) of test organism on non-selective media (e.g. Blood or Chocolate Agar)

QUALITY CONTROL

Frequency: Each shipment of each lot of reagent

Controls/Acceptable Results:

P. aeruginosa ATCC27853 - Oxidase Positive/purple
E. coli ATCC 25922 - Oxidase Negative/no color change

Corrective Actions for Unacceptable QC Results:

1. Repeat test using fresh subcultures of control organisms.
2. If still unsatisfactory, do not use for identification of isolates.
3. Arrange for replacement.

Documentation:

1. Record QC results on OXIDASE QC sheet.
2. Record unacceptable results and corrective actions on a QC DEVIATION form and submit to supervisor for review..
3. Supervisor will review QC monthly.

KCMC Biotechnology Laboratory, Microbiology	STANDARD OPERATING PROCEDURE	Effective Date 18 April 07	SOP-Number MIC.043
		Page 4 of 5	Date 18 April 2007
Title: OXIDASE TEST			

TEST PROCEDURE

1. Place a small piece of filter paper on a glass slide.
2. Add a drop of oxidase reagent to filter paper.
3. Using an inoculating loop or wooden applicator stick remove the colonies to be tested.
4. Smear the organism directly on to the reagent on the filter paper.

TEST READING

1. Examine the reaction area for development of purple color. Color will develop in 20 seconds.
2. Some organisms may produce a very light purple color or reaction may take longer than 20 sec. to develop (e.g. *Pasteurella* spp., mucoid *P. aeruginosa*).
3. A slight gray discoloration may develop with organisms that are oxidase negative – disregard.

INTERPRETATION

Positive – purple color development

Negative – no color change

RECORDING RESULTS

Enter results on POSITIVE CULTURE Form.

REFERENCE

Package Insert, BBL/BD DrySlide Oxidase, 2003. BD, Sparks, MD.

Kimberle, CC, Lauderdale, T. 2003. Reagents, Stains and Media. In: Manual of Clinical Microbiology, 8th Ed., ASM Press, Washington, DC.

APPENDIX A – Oxidase Test QC sheet

KCMC Biotechnology Laboratory, Microbiology	STANDARD OPERATING PROCEDURE	Effective Date 18 April 07	SOP-Number MIC.043
		Page 5 of 5	Date 18 April 2007
Title: OXIDASE TEST			

APPENDIX A

KCMC Biotechnology Laboratory
Microbiology

CONTROL ORGANISMS/ACCEPTABLE RESULTS:

Pseudomonas aeruginosa - Positive
E. coli - Negative

Quality Control
OXIDASE TEST

QC DATE	QC BY	LOT#	<i>P. aeruginosa</i> (+)	<i>E. coli</i> (-)	A/NA*

*A/NA = Acceptable/Not acceptable

DOCUMENT ALL CORRECTIVE ACTION ON QC

DEVIATION FORM

Supervisor Review :						

OXIDASE QC/QC SHEETS

Page 5 of 5